GB  
ES  
FR  
IT  
JP  
SC  
TC  
US

# USER INSTRUCTIONS

## **WIDEX STANDARD CHARGER**

WPT301 model

**WIDEX**

SOUND LIKE NO OTHER

|    |                                          |     |
|----|------------------------------------------|-----|
| GB | USER INSTRUCTIONS .....                  | 3   |
| ES | INSTRUCCIONES DE USO .....               | 24  |
| FR | MODE D'EMPLOI .....                      | 46  |
| IT | LIBRETTO D'ISTRUZIONI .....              | 68  |
| JP | 取扱説明書 .....                              | 89  |
| SC | 用户使用手册 .....                             | 111 |
| TC | 使用說明 .....                               | 131 |
| US | USER INSTRUCTIONS .....                  | 151 |
|    | FCC AND ISSED STATEMENTS .....           | 173 |
|    | COPYRIGHT WS AUDIOLOGY DENMARK A/S ..... | 176 |

# USER INSTRUCTIONS

|                                                  |           |
|--------------------------------------------------|-----------|
| <b>IMPORTANT SAFETY INFORMATION.....</b>         | <b>5</b>  |
| <b>YOUR CHARGER.....</b>                         | <b>7</b>  |
| Welcome.....                                     | 7         |
| The charger at a glance.....                     | 7         |
| Operating, storage and transport conditions..... | 8         |
| <b>OPERATING THE CHARGER.....</b>                | <b>9</b>  |
| Getting ready for use.....                       | 9         |
| Power supply.....                                | 9         |
| Power on.....                                    | 10        |
| Charging the hearing aid.....                    | 11        |
| Charging indicators.....                         | 12        |
| After charging.....                              | 14        |
| Overview of light indicators.....                | 15        |
| <b>CLEANING.....</b>                             | <b>17</b> |
| <b>TROUBLESHOOTING.....</b>                      | <b>18</b> |
| <b>REGULATORY INFORMATION.....</b>               | <b>20</b> |
| Intended use.....                                | 20        |
| EU directives.....                               | 20        |
| Directive 2014/53/EU.....                        | 20        |
| Information regarding disposal.....              | 20        |

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You can find technical data sheets and additional information on your charger on [www.widex.com](http://www.widex.com).

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If you are not comfortable with this, wait a few minutes for the hearing aid to reach room temperature before putting the hearing aid on.

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# CLEANING

The charger does not normally require cleaning. If necessary, disconnect the charger from the power supply and carefully wipe off the charger with a soft, dry cloth.

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To remove earwax from the charger, carefully wipe with a soft damp cloth.

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| Light indicators                                                                                                                                  | Possible cause                                         | Solution                                                                                                                                                                            |
|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Three or six indicators are flashing red every 2 seconds</p>  | <p>The charger cannot recognise the hearing aid(s)</p> | <p>Remove the incompatible hearing aid(s)</p>                                                                                                                                       |
|                                                                  | <p>The hearing aid has a critical error</p>            | <p>Remove the hearing aid, disconnect from the power supply for 10 seconds and then try again.<br/>If the charger is still flashing red, contact your hearing care professional</p> |
| <p>All indicators are flashing red in a sequence</p>             | <p>The charger has a critical error</p>                | <p>Remove the hearing aid, disconnect from the power supply for 10 seconds and then try again.<br/>If the charger is still flashing red, contact your hearing care professional</p> |

# REGULATORY INFORMATION

## Intended use

The WIDEX charger WPT301 model is to be used for charging Allure RIC and BTE hearing aids.

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Shipment of charger and hearing instruments should follow to local regulations.

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| Symbol                                                                           | Title/Description                                                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <b>Manufacturer</b><br>The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.                                                                                              |
|  | <b>Catalogue number</b><br>The product's catalogue (item) number.                                                                                                                                                                                                                 |
|  | <b>Consult instructions for use</b><br>The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.                                                                                                           |
|  | <b>Warning</b><br>Text marked with a warning symbol must be read before using the product.                                                                                                                                                                                        |
|  | <b>WEEE mark</b><br>“Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances. |

| Symbol                                                                            | Title/Description                                                                                                                                                              |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | <b>CE mark</b><br>The product is in conformity with the requirements set out in European CE marking directives.                                                                |
|  | <b>RCM mark</b><br>The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand markets. |
|  | <b>FCC logo</b><br>The product is compliant with relevant USA 47 CFR Rules.<br>FCC – Federal Communications Commission.                                                        |
|  | <b>Medical device</b><br>Indication that the device is a medical device.                                                                                                       |

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# USER INSTRUCTIONS

|                                                  |            |
|--------------------------------------------------|------------|
| <b>IMPORTANT SAFETY INFORMATION.....</b>         | <b>153</b> |
| <b>YOUR CHARGER.....</b>                         | <b>155</b> |
| Welcome.....                                     | 155        |
| The charger at a glance.....                     | 155        |
| Operating, storage and transport conditions..... | 156        |
| <b>OPERATING THE CHARGER.....</b>                | <b>157</b> |
| Getting ready for use.....                       | 157        |
| Power supply.....                                | 157        |
| Power on.....                                    | 158        |
| Charging the hearing aid.....                    | 159        |
| Charging indicators.....                         | 160        |
| After charging.....                              | 162        |
| Overview of light indicators.....                | 163        |
| <b>CLEANING.....</b>                             | <b>165</b> |
| <b>TROUBLESHOOTING.....</b>                      | <b>166</b> |
| <b>REGULATORY INFORMATION.....</b>               | <b>168</b> |
| Intended use.....                                | 168        |
| EU directives.....                               | 168        |
| Directive 2014/53/EU.....                        | 168        |
| Information regarding disposal.....              | 169        |

|                                  |            |
|----------------------------------|------------|
| Wireless radio transmission..... | 169        |
| <b>SYMBOLS.....</b>              | <b>171</b> |

# IMPORTANT SAFETY INFORMATION

Read these pages carefully before you start using your charger.

Keep the charger and other small parts away from children, pets, and mentally disabled people.

Never expose the charger to temperatures above 50 °C (122 °F), to prolonged sunlight, open flames or ignition sources.

Place the charger upright on a flat, stable surface when charging and with adequate space around the charger. Do not cover or place near flammable substances. This also applies when using a power bank.

Explosion hazard in explosive atmospheres! Do not use the charger in areas where there is a danger of explosions (e.g. mining area).

Do not drop. Internal components may get damaged.

Do not use if charger is broken.

Never try to open or repair the charger yourself. Contact your hearing care professional if you experience defects.

Never dry your device in the micro-wave or oven.

Only use a WPT301 charger with your hearing aids, if the WPT301 charger is compatible with your hearing aids. Use of a charger not recommended by WSAUD A/S may compromise the system or result in dangerous situations.

Never leave the charger in direct sunlight and do not immerse it in water.

Never insert other objects than your hearing aids in the charging wells. Do not insert incompatible connectors or foreign objects into the charger port.

The power adapter must comply with IEC 60601-1, IEC 62368-1 or equivalent safety standards acceptable in your region.

Always use the supplied USB cable with a power adapter with a USB-A type socket.

Use only a WSAUD charger with your hearing aids. Use of a charger not recommended by WSAUD A/S may compromise the system or result in dangerous situations.

# YOUR CHARGER

## Welcome

Thank you for choosing this product, which is intended to charge the supported hearing aid(s).

Read these instructions carefully before you start using this device.

### NOTE

Your charger may not look exactly like the one illustrated in this booklet. We reserve the right to make any changes we consider necessary.



## The charger at a glance

The charger has two charging wells and two sets of three light indicators that show the charging status of the hearing aid(s).

1. Charging wells

2. Light indicators

3. USB-C port

## Operating, storage and transport conditions

Your charger is designed for indoor use. It is sensitive to extreme conditions such as high temperatures, and it should not be operated, stored or transported in direct sunlight.

### Operate the charger within the following conditions:

|             | Minimum    | Maximum     |
|-------------|------------|-------------|
| Temperature | 0°C (32°F) | 30°C (86°F) |
| Humidity    | 10% rH     | 95% rH      |

### Store and transport the charger within the following conditions:

|                      | Minimum      | Maximum      |
|----------------------|--------------|--------------|
| Temperature          | -20°C (-4°F) | 50°C (122°F) |
| Humidity             | 10% rH       | 95% rH       |
| Atmospheric pressure | 750 mbar     | 1060 mbar    |


You can find technical data sheets and additional information on your charger on [www.wsa.com](http://www.wsa.com).

# OPERATING THE CHARGER



## Getting ready for use

### Power supply


You can connect your charger to different power supplies with 5 V DC, 500 mA output:

Use only the included USB cable.


Connect the cable to the charger:

1. USB-C port

## Power on


When the charger is connected to the power supply, the light indicators **will flash green in a sequence once** and then turn off:

The charger is now in standby mode and ready to charge your hearing aid.

## Charging the hearing aid

Gently place your hearing aid in the charging well. The hearing aid turns off automatically.

When connected to power, the charger can charge two hearing aids simultaneously; however, each hearing aid will be charged individually, depending on the state of the hearing aid.

This means that if a fully charged hearing aid is placed in the charger, it will not charge.

However, if the other hearing aid is *not* charged, it will charge.


Charging your hearing aid fully will take approximately 4 hours. If you charge your hearing aid for half an hour, you will have approximately 4 hours of use.

## Charging indicators

Once the charger has detected the hearing aid (s), the charging will begin. One light indicator will **glow green** to indicate that the hearing aid is low on battery and less than **33% charged**:


One light indicator will turn **solid green** and one will **glow green** to indicate that the hearing aid is more than **33% charged**:


Two light indicators will turn **solid green** and the third will **glow green** to indicate that the hearing aid is more than **67% charged**:

All three light indicators will **turn solid green** to indicate that the hearing aid is **100% charged**:

## **After charging**

Avoid removing the hearing aids from the charging wells by pulling the ear wire/tube, as this could damage the connection to the hearing aid.

Instead, put your fingers round the hearing aid when removing it from the charger.

When you remove the hearing aid from the charger, the hearing aid turns on after a few seconds.

Remember to turn the hearing aid and/or the charger off if you are not going to use them right away.

Immediately after charging, the hearing aid may temporarily feel warm.

If you are not comfortable with this, wait a few minutes, for the hearing aid to reach room temperature, before putting the hearing aid on.

## Overview of light indicators

### Light indicators

### Status

All indicators will flash green in a sequence once and then turn off

The charger is connected to power.

One single indicator will glow green

The hearing aid (s) is charging.  
Low battery, less than 33%

One indicator will turn solid green and one will glow green

The hearing aid (s) is charging.  
Battery more than 33%

Two indicators will turn solid green and the third will glow green

The hearing aid (s) is charging.  
Battery more than 67%

All indicators will turn solid green

The hearing aid (s) is fully charged.  
Battery 100% charged

## Light indicators

## Status

Indicators will flash red three times every 2 seconds

A hearing aid error occurred during charging and charging has stopped.  
See Troubleshooting

All indicators will flash red in a sequence

A charger error occurred during charging and charging has stopped.  
See Troubleshooting

All indicators will flash amber in a sequence

The temperature of the hearing aid (s) has become too high during charging and charging has been paused.  
See Troubleshooting

# CLEANING

The charger does not generally require cleaning. If necessary, disconnect the charger from the power supply and carefully wipe off the charger with a soft dry cloth

If you notice any dust in the charging well, gently blow it clean. If this is not enough, gently clean the charging well with a cotton bud.

To remove earwax from the charger, carefully wipe with a soft cloth with a bit of mild household detergent.

Never use sharp/pointy objects to clean the charger.

# TROUBLESHOOTING

These pages contain advice on what to do if your charger does not work as expected.

| Light indicators                                                                                                                             | Possible cause                                                                                                            | Solution                                                                                                                                                                     |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No light indicators                                                                                                                          | The charger is not connected                                                                                              | Make sure that the USB cable is connected to a power adapter and then plugged into a wall socket                                                                             |
|                                                                                                                                              | The hearing aid is not placed correctly in the charger                                                                    | Place the hearing aid correctly in the charger                                                                                                                               |
| Three or six indicators are flashing amber in a sequence<br> | The temperature of your surroundings or of the hearing aid has exceeded a critical limit and the charging has been paused | The charging will continue when the temperature of your surroundings or of the hearing aid has dropped.<br>If this continues, find a cooler place to charge your hearing aid |

| Light indicators                                                                                                                                  | Possible cause                                         | Solution                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Three or six indicators are flashing red every 2 seconds</p>  | <p>The charger cannot recognise the hearing aid(s)</p> | <p>Remove the incompatible hearing aid(s)</p>                                                                                                                                   |
|                                                                  | <p>The hearing aid has a critical error</p>            | <p>Remove the hearing aid, disconnect from the power supply for 10 seconds and then try again. If the charger is still flashing red, contact your hearing care professional</p> |
| <p>All indicators are flashing red in a sequence</p>             | <p>The charger has a critical error</p>                | <p>Remove the hearing aid, disconnect from the power supply for 10 seconds and then try again. If the charger is still flashing red, contact your hearing care professional</p> |

# REGULATORY INFORMATION

## Intended use

The Widex charger WPT301 model is to be used for charging Allure RIC and BTE hearing aids.

## EU directives

### Directive 2014/53/EU

Hereby, WSAUD A/S declares that the WPT301 is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

The WPT301 contains a radio transmitter operating at: 6.78 MHz, 13.4 dB $\mu$ A/m @ 3 m.

A copy of the Declaration of Conformity according to 2014/53/EU can be found at

[www.widex.com/doc](http://www.widex.com/doc)

N26346

## Information regarding disposal

Do not dispose of hearing aids, hearing aid accessories and chargers with ordinary household waste.

Hearing aids, hearing aid accessories and chargers should be disposed of at sites intended for electrical waste and electronic equipment, or given to your hearing care professional for safe disposal.

For returning or shipping the devices, please contact your hearing care professional.

## Wireless radio transmission

The product complies with the following EMC/EMI standards:

| Standard                  | Test type                                            | Note                                                         |
|---------------------------|------------------------------------------------------|--------------------------------------------------------------|
| 47 CFR Part 15, Subpart B | RF emissions                                         | US FCC requirements for unintentional radiators              |
| IEC 60601-1-2             | General EMC including RF emissions, RF immunity, ESD | Medical Electrical Equipment, Electro-magnetic compatibility |

| Standard     | Test type   | Note                                                                                                                 |
|--------------|-------------|----------------------------------------------------------------------------------------------------------------------|
| IEC 60118-13 | RF immunity | International standard for hearing instruments to ensure adequate immunity radio interference from mobile telephones |

# SYMBOLS

Symbols commonly used by WSAUD A/S in medical device labelling (labels/IFU/etc.)

| Symbol                                                                            | Title/Description                                                                                                                                                                                                                                                                 |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <b>Manufacturer</b><br>The product is produced by the manufacturer whose name and address are stated next to the symbol. If appropriate, the date of manufacture may also be stated.                                                                                              |
|  | <b>Catalogue number</b><br>The product's catalogue (item) number.                                                                                                                                                                                                                 |
|  | <b>Consult instructions for use</b><br>The user instructions contain important cautionary information (warnings/precautions) and must be read before using the product.                                                                                                           |
|  | <b>Warning</b><br>Text marked with a warning symbol must be read before using the product.                                                                                                                                                                                        |
|  | <b>WEEE mark</b><br>“Not for general waste”. When a product is to be discarded, it must be sent to a designated collection point for recycling and recovering to prevent the risk of harm to the environment or human health as a result of the presence of hazardous substances. |

| Symbol                                                                           | Title/Description                                                                                                                                                              |
|----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | <b>CE mark</b><br>The product is in conformity with the requirements set out in European CE marking directives.                                                                |
|  | <b>RCM mark</b><br>The product complies with electrical safety, EMC and radio spectrum regulatory requirements for products supplied to the Australian or New Zealand markets. |
|  | <b>FCC logo</b><br>The product is compliant with relevant USA 47 CFR Rules.<br>FCC – Federal Communications Commission.                                                        |
|  | <b>Medical device</b><br>Indication that the device is a medical device.                                                                                                       |

# FCC AND ISED STATEMENTS

## Federal Communications Commission Statement

This device complies with part 15/part 18 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

### NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**NOTE:**

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Changes or modifications to the equipment not expressly approved by WSAUD A/S could void the user's authority to operate the equipment.

**ISED Statement / Déclaration d'ISED**

CAN RSS-216/CNR-216

CAN ICES-1/NMB-1

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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#### ISED RADIATION EXPOSURE STATEMENT:

This equipment complies with ISED RF radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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# **COPYRIGHT WS AUDIOLOGY DENMARK A/S**

Copyright © WS Audiology Denmark A/S. All rights reserved. Any unauthorized redistribution or use of the software in source or binary form is prohibited and will be prosecuted when discovered.

**WSAUD A/S**

Nymoellevej 6, DK-3540 Lyngby, Denmark  
[www.widex.com](http://www.widex.com)

**Manual no.:**  
9 514 0816 094 01  
**Issue:**  
2023-09